

VAHTS Maxi Unique Dual Index DNA Adapters Set 2 for Illumina

N34202



Nanjing Vazyme Biotech Co.,Ltd.

Tel: +86 25-83772625

Email: info.biotech@vazyme.com

Web: www.vazyme.com

Loc: Red Maple Hi-tech Industry Park, Nanjing, PRC

Follow Us



Instruction for Use

Version 22.1

Contents

01/Product Description	02
02/Components	02
03/Storage	02
04/Applications	02
05/Notes	02
06/Library Structure and Sequences	02
07/Appendix	05

01/Product Description

VAHTS Maxi Unique Dual Index DNA Adapters Set 2 for Illumina is a dedicated kit for multi-sample DNA library preparation for sequencing using Illumina high-throughput sequencing platforms. The kit contains 96 kinds of VAHTS Unique Dual Index (UDI) DNA Adapters, which enables dual verification through two completely independent indexes and minimizes index hopping and misassignment, thus ensuring that the reads in the final analysis represent the samples to the maximum extent.

All the Adapters provided in the kit have undergone rigorous quality control and functional testing to ensure the optimal stability and repeatability of library preparation.

02/Components

Components	N34202-01 (384 rxns)
VAHTS Maxi Unique Dual Index DNA Adapter for Illumina (UDIA 097 - UDIA 0192)	20 µl each

03/Storage

Store at -30 ~ -15°C and transport at ≤0°C.

04/Applications

VAHTS Maxi Unique Dual Index DNA Adapters Set 2 for Illumina is a dedicated kit for Vazyme ND series library preparation kit. It is compatible with various template type: genomic DNA, cfDNA, ctDNA and FFPE DNA. It is suitable for multi-sample Dual indexed DNA library preparation and can effectively prevent crosstalk between samples.

05/Notes

For research use only. Not for use in diagnostic procedures.

1. The use amount of DNA Adapter depends on the Input DNA amount. Refer to Appendix for specific concentration.
2. Please do not premix Adapter, Ligation Buffer and Ligase, or it may result in Adapter self-ligation.
3. Please do not keep this product at room temperature, or the ligation efficiency may decrease.

06/Library Structure and Sequences

The structure of a DNA library prepared by using VAHTS Maxi Unique Dual Index DNA Adapters for Illumina is as follows:

5' - **DNA Adapter X** - Insert DNA Sequence - **DNA Adapter X** - 3'

Index sequences are as follows:

Adapter name	i5 index	i5 index	i7 index
	(HiSeq 2000/2500, MiSeq, NovaSeq v1.0)	(HiSeq 3000/4000, NextSeq, MiniSeq, NovaSeq v1.5)	(all Illumina systems)
UDIA-097	CGACCATT	AATGGTCG	GCCTATCA
UDIA-098	GATAGCGA	TCGCTATC	CTTGGATG
UDIA-099	AATGGACG	CGTCCATT	AGTCTCAC
UDIA-100	CGCTAGTA	TACTAGCG	CTCATCAG
UDIA-101	TCTCTAGG	CCTAGAGA	TGTACCGT
UDIA-102	ACATTGCG	CGCAATGT	AAGTCGAG
UDIA-103	TGAGGTGT	ACACCTCA	CACGTTGT
UDIA-104	AATGCCTC	GAGGCATT	TCACAGCA
UDIA-105	CTGGAGTA	TACTCCAG	CTACTTGG
UDIA-106	GTATGCTG	CAGCATAC	CCTCAGTT
UDIA-107	TGGAGAGT	ACTCTCCA	TCCTACCT
UDIA-108	CGATAGAG	CTCTATCG	ATGGCGAA
UDIA-109	CTCATTGC	GCAATGAG	CTTACCTG
UDIA-110	ACCAGCTT	AAGCTGGT	CTCGATAC
UDIA-111	GAATCGTG	CACGATTC	TCCGTGAA
UDIA-112	AGGCTTCT	AGAAGCCT	TAGAGCTC
UDIA-113	CAGTTCTG	CAGAACTG	TGACTGAC
UDIA-114	TTGGTGAG	CTCACCAA	TAGACGTG
UDIA-115	CATTGCGT	ACCGAATG	CCGGAATT
UDIA-116	TGTGAAGC	GCTTCACA	CTCCTAGA
UDIA-117	TAAGTGGC	GCCACTTA	CAACGGAT
UDIA-118	ACGTGATG	CATCACGT	TGGCTATC
UDIA-119	GTAGAGCA	TGCTCTAC	CGGTCATA
UDIA-120	GTCAGTTG	CAACTGAC	TCCAATCG
UDIA-121	ATTCGAGG	CCTCGAAT	GAGCTTGT
UDIA-122	GATACTGG	CCAGTATC	GAAGGTTT
UDIA-123	GCCTTGTT	AACAAGGC	ATCTCGCT
UDIA-124	TTGGTCTC	GAGACCAA	AGTTACGG
UDIA-125	CCGACTAT	ATAGTCGG	GTGTCTGA
UDIA-126	GTCTTAAG	CTTAGGAC	TGACTTCG
UDIA-127	ACCAATGC	GCATTGGT	TGGATCAC
UDIA-128	GATGCACT	AGTGCATC	ACACCACT
UDIA-129	GCTGGATT	AATCCAGC	CAGGTTAG
UDIA-130	ATGGTTGC	GCAACCAT	AGTTGGCT
UDIA-131	CAGAATCG	CGATTCTG	TCAACTGG
UDIA-132	GAACGCTT	AAGCGTTC	CTGCAGTT
UDIA-133	TCGAACCA	TGTTTCGA	ACACGGTT
UDIA-134	CTATCGCA	TGCGATAG	AATACGCG
UDIA-135	TACGGTTG	CAACCGTA	TGCGAACT
UDIA-136	GAGATGTC	GACATCTC	GCTGACTA
UDIA-137	CTTACAGC	GCTGTAAG	GTGGTGTT
UDIA-138	AGGAGGAA	TTCTCTCT	GTGCTTAC
UDIA-139	GACGAATG	CATTGCTC	TCAAGGAC
UDIA-140	GAAGAGGT	ACCTCTTC	TGAACCTG

N34202

UDIA-141	CGTCAATG	CATTGACG	AGTGTGGG
UDIA-142	TACCAGGA	TCCTGGTA	GTACTCTC
UDIA-143	CGTACGAA	TTCGTACG	CCGTATCT
UDIA-144	GACTTAGG	CCTAAGTC	CGAAGAAC
UDIA-145	AGTGCAGT	ACTGCACT	AGCGGAAT
UDIA-146	TTGATCCG	CGGATCAA	GTGAGCTT
UDIA-147	TGCCATTG	GAATGGCA	CGTGATCA
UDIA-148	CTTGCTGT	ACAGCAAG	TCGCATTG
UDIA-149	CCTACTGA	TCAGTAGG	TGACGCAT
UDIA-150	CCAAGTTG	CAACTTGG	CCGATGTA
UDIA-151	TGATCGGA	TCCGATCA	TTCCGAGT
UDIA-152	TAGTTGCG	CGCAACTA	ACGACAGA
UDIA-153	GTCTGATC	GATCAGAC	AGCTTGAG
UDIA-154	CGTTATGC	GCATAACG	GAGTGGTT
UDIA-155	GCTCTGTA	TACAGAGC	GCTGTAAG
UDIA-156	TTACCGAG	CTCGGTAA	CCAAGACT
UDIA-157	GCCATAAC	GTTATGGC	ATTGCGTG
UDIA-158	CTCAGAGT	ACTCTGAG	CTGAAGCT
UDIA-159	CGAGACTA	TAGTCTCG	TAACGAGG
UDIA-160	TGTGCGTT	AACGCACA	TCGTCTCA
UDIA-161	TTCAGGAG	CTCCTGAA	TTCTGTGT
UDIA-162	GACTATGC	GCATAGTC	CGTTGAGT
UDIA-163	AGGTTCTG	TCGAACCT	AGTCGCTT
UDIA-164	AGTCTGTG	CACAGACT	TAGGTAGG
UDIA-165	ACCTAAGG	CCTTAGGT	CAGGAGAT
UDIA-166	TGCAGGTA	TACCTGCA	CATCGTGA
UDIA-167	AAGGACAC	GTGTCCTT	TGTTGTGG
UDIA-168	CAACCTAG	CTAGGTTG	ACAGACCT
UDIA-169	CTGACACA	TGTGTCAG	GTCTTTCT
UDIA-170	ACTCGTTG	CAACGAGT	TGATACGC
UDIA-171	AGTCCTTA	TAGGAGCT	CTGTGTTG
UDIA-172	TACATCGG	CCGATGTA	AACGTGGA
UDIA-173	CACAAGTC	GACTTGTG	GTGCGCAT
UDIA-174	CGGATTGA	TCAATCCG	AACGACGT
UDIA-175	AGTCGACA	TGTCGACT	CGTATTCT
UDIA-176	GTCTCCTT	AAGGAGAC	AGCAAGCA
UDIA-177	GAGATACG	CGTATCTC	TGTTCCAG
UDIA-178	ATCGGTGT	ACACCGAT	CTCCATGT
UDIA-179	TTCGCAA	TTGCGAGA	CGCTTTGT
UDIA-180	TCTAACGC	GCGTTAGA	ATAAGGCG
UDIA-181	CAATCGAC	GTCGATTG	TGTCTGCT
UDIA-182	GAGGACTT	AAGTCCTC	CGCTTAAC
UDIA-183	TGGAGTTG	CAACTCCA	GATCCATG
UDIA-184	CTAGGCAT	ATGCCTAG	ACCTCTGT
UDIA-185	CTCTACTC	GAGTAGAG	GCCACTTA
UDIA-186	AGAAGCGT	ACGCTTCT	ACCTGACT
UDIA-187	TCGAAGGT	ACCTTCGA	GTAAAGGC
UDIA-188	GTGCGTAA	TTACCGAC	ATGCCAAC
UDIA-189	ACGATGAC	GTCATCGT	AGAGGTTG
UDIA-190	TCCGTATG	CATACGGA	ACCATCCA
UDIA-191	CTAGGTGA	TCACCTAG	GTGGATAG
UDIA-192	CATTGCCT	AGGCAATG	CTGAGATC

07/Appendix

Table 1. Recommended Adapter Dilution Ratios and Volumes for Library Preparation Featuring Mechanical Fragmentation

Input DNA	UDI Adapter Dilution Ratio	Usage Volume
1 ng	1:90	5 µl
10 ng	1:15	5 µl
50 ng	1:4	5 µl
100 ng	1:2	5 µl
≥500 ng	No dilution	5 µl

▲ Library preparation kit: VAHTS Universal DNA Library Prep Kit for Illumina V3 (Vazyme #ND607).

Table 2. Recommended Adapter Dilution Ratios and Volumes for Library Preparation Featuring Enzymatic Fragmentation

Input DNA	UDI Adapter Dilution Ratio	Usage Volume
1 ng	1:70	5 µl
100 ng	1:2	5 µl
1 µg	No dilution	5 µl

▲ Library preparation kit: VAHTS Universal Plus DNA Library Prep Kit for Illumina (Vazyme #ND617).

Table 3. Recommended Adapter Dilution Ratios and Volumes for Transcript Library Preparation

Input RNA	UDI Adapter Dilution Ratio	Usage Volume
10 ng	1:10	0.5 µl
100 ng	1:10	1 µl
1 µg	1:10	3.5 µl

▲ Library preparation kit: VAHTS Universal V8 RNA-seq Library Prep Kit for Illumina (Vazyme #NR605).