AccurSTART Reverse Transcriptase (Glycerol-free)



Version 23.1



Product Description

AccurSTART Reverse Transcriptase (Glycerol-free) is a lyophilizable version of AccurSTART Reverse Transcriptase (Vazyme #RW201). It can be applied to downstream lyophilization technology while maintaining the great reverse transcription performance and stability of Vazyme #RW201. This product does not contain excipients, please add your own as needed.

Components

| Components | RWL201-01 (10,000 U) | RWL201-02 (40,000 U) |
|---|-------------------------|-------------------------|
| 5 × AccurSTART Buffer | 200 µl | 800 µl |
| AccurSTART Reverse Transcriptase (Glycerol-free) (200 U/μl) | 50 μl | 200 μΙ |

Storage

Store at -30 ~ -15°C and transport at ≤0°C.

Applications

It is applicable for one-step RT-qPCR reactions.

Source

A recombinant E. coli strain carrying modified M-MLV (H-) reverse transcriptase gene.

Unit Definition

One unit (U) is defined as the amount of enzyme that incorporates 1 nmol of dTTP into acid-insoluble material in 10 min at 37°C, with Poly(rA)·Oligo (dT) as the template/primer.

Notes

For research use only. Not for use in diagnostic procedures.

- 1. Please keep the experiment area clean; Wear disposable gloves and masks; Use RNase-free consumables such as centrifuge tubes and pipette tips.
- 2. Keep RNA on ice to avoid degradation.
- 3. High quality RNA templates are recommended to achieve high efficiency reverse transcription.