

Ribo-clean rRNA Depletion Kit Mega (Prevalent Species)

RN416



Instruction for Use

Version 24.1

Contents

01/Product Description	02
02/Target Species	02
03/Components	02
04/Storage	02
05/Applications	02
06/Self-prepared Materials	03
07/Notes	03
07-1/About Storage	03
07-2/About RNA Sample Preparation	03
07-3/About RNA Magnetic Beads	04
07-4/About Operation	04
08/Mechanism & Workflow	05
09/Experiment Process	06
09-1/Probe Hybridization with RNA	06
09-2/RNase H Digestion	06
09-3/DNase I Digestion	07
09-4/Ribosomal-depleted RNA Purification	07
10/FAQ & Troubleshooting	08

For Research Use Only. Not for use in diagnostic procedures.

01/Product Description








Ribo-clean rRNA Depletion Kit Mega (Prevalent Species) is designed for rRNA depletion with 0.01 - 1 µg total RNA from human, rat, mouse, and macaca. The rRNA (including 5S, 5.8S, 18S, 28S rRNA, mt12S, mt16S rRNA, and 45S rRNA) is depleted after rRNA and probe hybridization, RNase H digestion, and DNase I digestion, while mRNA and other non-coding RNA (e.g., lncRNA) are retained for later analysis. The kit is applicable to intact or partially degraded RNA samples, and the optimized reaction system improves the depletion efficiency and species compatibility.

02/Target Species

Cat. No.	Target Species*
RN416	Human/Rat/Mouse/ <i>Cricetulus griseus</i> / <i>Macaca mulatta</i> / <i>Macaca fascicularis</i>

* Species coverage can be found at <http://cloud.vazyme.com:83/application-tool>.

03/Components

Components	RN416-01 (6 rxns)	RN416-02 (24 rxns)	RN416-03 (96 rxns)
 rRNA Probe Mega (Prevalent Species)	18 µl	72 µl	288 µl
 Probe Buffer	18 µl	72 µl	288 µl
 RNase H Buffer 2	18 µl	72 µl	288 µl
 RNase H Mix	12 µl	48 µl	192 µl
 DNase I Buffer	174 µl	696 µl	4 × 696 µl
 DNase I	6 µl	24 µl	96 µl
 Nuclease-free ddH ₂ O	1 ml	1 ml	4 × 1 ml

▲ The colors in the table correspond to the cap color of each component.

04/Storage

Store at -30 ~ -15°C and ship at ≤0°C.

05/Applications

Ribo-clean rRNA Depletion Kit Mega (Prevalent Species) is designed for rRNA depletion (including 5S, 5.8S, 18S, 28S rRNA, mt12S, mt16S rRNA, 45S rRNA) with 0.01 - 1 µg total RNA from human, rat, mouse, and macaca. It is compatible with partially degraded RNA samples, and the depleted product is applicable to RNA library preparation and other experiments. The amount of mRNA varies significantly across different samples, and the initial total RNA inputs can be adjusted according to downstream applications. For RNA library preparation, VAHTS Universal V8 RNA-seq Library Prep Kit for Illumina (Vazyme #NR605), VAHTS Universal V10 RNA-seq Library Prep Kit for Illumina (Vazyme #NR606), VAHTS Universal V10 RNA-seq Library Prep Kit (Premixed Version) (Vazyme #NR616) are recommended.

06/Self-prepared Materials

- ◇ RNA Evaluation
Agilent RNA 6000 Pico Kit (Agilent #5067-1513).
- ◇ RNA Purification
VAHTS RNA Clean Beads (Vazyme #N412).
- ◇ Library Preparation
VAHTS Universal V8 RNA-seq Library Prep Kit for Illumina (Vazyme #NR605),
VAHTS Universal V10 RNA-seq Library Prep Kit for Illumina (Vazyme #NR606),
VAHTS Universal V10 RNA-seq Library Prep Kit (Premixed Version) (Vazyme #NR616).
- ◇ Other Materials
80% ethanol (freshly prepared with Nuclease-free ddH₂O), Nuclease-free ddH₂O;
Nuclease-free PCR tubes; Agilent 2100 Bioanalyzer or other equivalent products, PCR
instrument, magnetic rack (Vazyme #CM101), etc.

07/Notes

07-1/About Storage

1. The kit contains multiple enzymes, which must be stored at -30 ~ -15°C. When in use, it should be placed on ice and stored as required after use, or the enzyme activity may be reduced.
2. To avoid repeated freeze-thaw cycles, please aliquot and store the remaining reagent after first use.

07-2/About RNA Sample Preparation

1. To ensure depletion efficiency, RNA samples should be free of salts (e.g., Mg²⁺ or guanidine salts) or organics (e.g., phenol or ethanol), or another purification is required.
2. To avoid DNA contamination, RNA samples can be treated with DNase I to remove DNA.
3. Dilute RNA with Nuclease-free ddH₂O to 9 µl and avoid prolonged exposure on ice to avoid RNA degradation.
4. If a low RNA concentration results in an initial volume >9 µl, the RNA can be concentrated by lyophilization, ethanol precipitation, column purification or magnetic bead purification (VAHTS RNA Clean Beads, Vazyme #N412).
5. For RNA-Seq applications, the total RNA input >100 ng is recommended to increase library complexity.

07-3/About RNA Magnetic Beads


1. Equilibrate magnetic beads to room temperature and mix well before use, or the recovery efficiency may be affected.
2. Transfer the supernatant when the solution becomes clear and do not disturb the beads.
3. Use 80% ethanol (freshly prepared with Nuclease-free ddH₂O) to rinse magnetic beads for ribosomal-depleted RNA purification; otherwise it may cause RNA loss and library preparation failure. Please discard as much of the supernatant as possible to minimize impurity residues after second rinse.
4. Before elution, ensure that the magnetic beads are sufficiently dry (with the surface changing from shiny brown to matte brown) to prevent ethanol residue from affecting downstream experiments. However, it is important to avoid over-drying as it can reduce the RNA recovery efficiency.

07-4/About Operation

1. It is recommended to use pipette tips with filters and change the pipette tips for different samples.
2. Be sure to wear gloves during operation, and please change gloves after touching equipment or other work areas outside the RNase-free zone.
3. All reagents must be capped immediately after use to avoid contamination.
4. Briefly centrifuge enzyme components before use to avoid adhesion to the tube walls and caps.

08/Mechanism & Workflow

1. rRNA probe hybridization



2. RNase H digestion



3. DNase I digestion



— rRNA
— mRNA or nc-RNA
— probe

Fig 1. Schematic Diagram of rRNA Depletion




Fig 2. Workflow of rRNA Depletion


09/Experiment Process

09-1/Probe Hybridization with RNA

- 1. Dilute 0.01 - 1 µg total RNA to a final volume of 9 µl with Nuclease-free ddH₂O in the Nuclease-free PCR tube, and keep it on ice for later use.

▲ Take out the required components from -30 ~ -15°C in advance for next step and place on ice for later use.

- 2. Prepare the following reaction mix in a Nuclease-free PCR tube:

Components	Volume
rRNA Probe Mega (Prevalent Species)	3 µl 
Probe Buffer	3 µl 
Total RNA	9 µl
Total	15 µl

- 3. Mix well by pipetting 10 times, and briefly centrifuge to collect the reaction mix at the bottom of the tube.

▲ When multiple samples are processed simultaneously, the rRNA Probe Mega (Prevalent Species) and Probe Buffer can be pre-mixed, and then aliquoted into each PCR tube. It is recommended to prepare 1.1 times the actual number of reactions to compensate for losses.



- 4. Place the PCR tube in the PCR instrument and perform the following program:

Temperature	Time
95°C	2 min
95 ~ 37°C	0.1°C/sec
37°C	5 min
4°C	Hold

- 5. Remove the tube immediately when the reaction is completed. Briefly centrifuge to collect the reaction mix at the bottom of the tube, and place it on ice for the next step.

09-2/RNase H Digestion

- 1. Prepare the following reaction mix on ice:

Components	Volume
RNase H Buffer 2	3 µl 
RNase H Mix	2 µl 
Previous step product	15 µl
Total	20 µl

- 2. Mix well by pipetting 10 times, and briefly centrifuge to collect the reaction mix at the bottom of the tube.



- 3. Place the PCR tube in the PCR instrument and perform the following program:

Temperature	Time
50°C	15 min
4°C	Hold

- Briefly centrifuge to collect the reaction mix at the bottom of the tube, and place it on ice for the next step.

09-3/DNase I Digestion

- Prepare the following reaction mix on ice:

Components	Volume
DNase I Buffer	29 μ l 
DNase I	1 μ l 
Previous step product	20 μ l
Total	50 μ l

- Mix well by pipetting 10 times, and briefly centrifuge to collect the reaction mix at the bottom of the tube.
- Place the PCR tube in the PCR instrument and perform the following program:

Temperature	Time
37°C	10 min
4°C	Hold

- Briefly centrifuge to collect the reaction mix at the bottom of the tube, and place it on ice for the next step.

09-4/Ribosomal-depleted RNA Purification

- Mix VAHTS RNA Clean Beads (Vazyme #N412) well by vortexing. Add 110 μ l (2.2 \times) to the RNA sample from the previous step, and mix well by pipetting 10 times.
- Incubate on ice for 15 min to allow RNA to bind to the magnetic beads.
- Place the tube on the magnetic rack, and discard the supernatant when the solution becomes clear (about 5 min).
- Keep the tube on the magnetic rack, and wash the beads with 200 μ l of 80% ethanol (freshly prepared with Nuclease-free ddH₂O). Incubate at room temperature for 30 sec, then carefully discard the supernatant.
- Repeat step 4 once.
- Keep the tube on the magnetic rack, and air-dry the beads for 5 - 10 min.

▲ Do not disturb the magnetic beads when adding 80% ethanol.

▲ Use a 10 μ l pipette to discard the residual supernatant.

▲ Avoid over-drying of the magnetic beads, which may reduce recovery efficiency.

7. RNA elution:

Option A - Reverse transcription application:

Remove the tube from the magnetic rack, add 20 µl of Nuclease-free ddH₂O, and mix well. Incubate at room temperature for 2 min. Place the tube on the magnetic rack, and carefully transfer 18 µl of the supernatant to a new Nuclease-free PCR tube when the solution becomes clear (about 5 min). Store at -85 ~ -65°C for later use.

Option B - RNA library preparation application:

Remove the tube from the magnetic rack, add 18 µl of Frag/Prime Buffer 2 from VAHTS Universal V10 RNA-seq Library Prep Kit for Illumina (Vazyme #NR606), and mix well. Incubate at room temperature for 2 min. Place the tube on the magnetic rack, and carefully transfer 16 µl of the supernatant to a new Nuclease-free PCR tube when the solution becomes clear (about 5 min), and immediately proceed to library preparation.

10/FAQ & Troubleshooting

◇ If rRNA depletion is required for multiple species, can multiple probes be used simultaneously?

Multiple probes can be used simultaneously for rRNA depletion across various species. Follow the recommended amount of probe in the instruction manual and reduce the volume of RNA input to maintain a total volume of 15 µl.

Take multi-species rRNA depletion (Vazyme #RN416+RN417) for example:

Components	Volume
rRNA Probe Mega (Prevalent Species)	3 µl
rRNA Probe Mega (Bacteria)	3 µl
Probe Buffer	3 µl
Total RNA	6 µl
Total	15 µl

▲ Only the reaction system in 09-1/Probe Hybridization with RNA is changed, while all other operations remain unchanged

▲ If the initial amount of RNA is insufficient due to reduced volume, the RNA can be concentrated by VAHTS RNA Clean Beads (Vazyme #N412).

◇ How to know the species coverage and compatible products?

Please refer to <http://cloud.vazyme.com:83/application-tool> to get the species coverage. The species coverage is not the same as the depletion efficiency and is for reference only.

◇ How to store the purified product?

The purified product is prone to degradation due to low concentration. Proceed to downstream experiments as soon as possible, or store at -85 ~ -65°C.

◇ If the purified product is used for library preparation, but it is eluted with Nuclease-free ddH₂O, what's the operation?

Purify the RNA again with VAHTS RNA Clean Beads (Vazyme #N412), and finally elute with Frag/Prime Buffer.

◇ If the library concentration is too low, how to improve it?

Library concentration is related to total RNA quality including the integrity and input amounts. If it is not possible to extract qualified RNA samples, the following compensatory methods may be taken:

1. Increase total RNA inputs: up to 1 µg.
2. Prepare several duplicate samples and merge them after purification.



Vazyme Biotech Co.,Ltd.

www.vazyme.com

400-600-9335 (China) +86 400-168-5000 (Global)

support@vazyme.com