

# Ultra-Universal TOPO Cloning Kit

C603

Version 23.1



## Product Description

Ultra-Universal TOPO Cloning Kit is the fourth-generation TOPO cloning kit. The kit uses upgraded Topoisomerase with the optimal buffer and has higher enzyme activity. The addition of the second-generation blunting factor and vector self-ligation inhibitor makes this product compatible with both TA cloning and blunt-end cloning, and increases the positive rate of cloning.

## Components

Components	C603-01 (20 rxns)	C603-02 (40 rxns)
■ 5 × Ultra-Universal TOPO Cloning Mix <sup>a</sup>	40 μl	2 × 40 μl
■ 500 bp Control insert (20 ng/μl)	10 μl	20 μl
■ M13 Primer Mix (10 μM) <sup>b</sup>	200 μl	400 μl

a. It contains Topoisomerase and pCE3 Blunt Vector (Amp resistance vector).

b. It contains M13 Forward Primer and M13 Reverse Primer.

## Storage

Store at -30 ~ -15°C and transport at ≤0°C.


## Applications

It is applicable for the ligation of blunt-end or A-tailed PCR products.

## Notes

For research use only. Not for use in diagnostic procedures.

## Mechanism & Workflow



Workflow of Ultra-Universal TOPO Cloning Kit

Figure A: TA Cloning

- Add A-tailed PCR products of *Taq* DNA polymerase (such as Vazyme's *Taq* series) into 5 × Ultra-Universal TOPO Cloning Mix, then react at 20 ~ 37°C for 5 min.
- Blunting factor in Mix removes the adenine at the end of the amplification product to form blunt-end products.
- The 5'-OH of blunt-end products attack the phosphate bond between the topoisomerase and the vector. Then, the topoisomerase is released, and the vector forms a circular recombinant with the blunt-end product.

Figure B: Blunt-end Cloning

- Add blunt-end PCR products of high-fidelity DNA polymerase (such as Vazyme's *Phanta* series) into 5 × Ultra-Universal TOPO Cloning Mix, then react at 20 ~ 37°C for 5 min.
- The 5'-OH of blunt-end products attack the phosphate bond between the topoisomerase and the vector. Then, the topoisomerase is released, and the vector forms a circular recombinant with the blunt-end product.

## Experiment Process

### 1. Preparation of PCR products

- Primer requirements: The 5' end of the primers cannot be phosphorylated.
- Enzyme selection: It is recommended to use Vazyme's *Phanta* or *Taq* series of DNA polymerase products.



c. Product requirements: Please ensure the integrity of PCR amplification products; after amplification reaction, it is recommended to perform electrophoresis to detect the yield and quality of products, and then perform gel extraction and purification.

## 2. Ligation Reaction

### a. Reaction system

**The optimal mass of insert required = [0.03 × number of base pairs] ng**

For example, the length of insert is 1,000 bp: the optimal mass of insert required = [0.03 × 1,000] ng = 30 ng. This product has a broad range of input compatibility for inserts. The recommended amount is as follows:

Inserts Size	Recommended Mass
0.05 - 1 kb	5 - 30 ng
1 - 2 kb	30 - 60 ng
2 - 5 kb	60 - 150 ng
>5 kb	>150 ng

▲ It is recommended to use Nanodrop, Onedrop, etc. for concentration measurement.

Components	Volume
5 × Ultra-Universal TOPO Cloning Mix	2 µl
Purified PCR product	1 - 8 µl
ddH <sub>2</sub> O	to 10 µl

- b. Mix by flicking the bottom of the tube, then centrifuge briefly to collect all the liquid to the bottom of the centrifuge tube.
- c. This kit has high compatibility with reaction temperature, it only needs to react at 20 ~ 37°C for 5 min. It is recommended to react at 25°C for 5 min, and use the PCR machine to control the temperature.
- d. After the reaction, place the centrifuge tube on ice.

## 3. Transformation

Add 10 µl of the ligation product to 100 µl of conventional chemically competent cells, and refer to the competent instructions for transformation conditions.

## 4. Positive Clone Identification

**a. Colony/Bacteria Liquid PCR Identification:** Pick a single clone into 10 µl ddH<sub>2</sub>O and mix well as template; 2 × Rapid Taq Master Mix (Vazyme #P222) is recommended.

### Reaction System

Components	Volume
2 × Rapid Taq Master Mix	10 µl
M13 Primer Mix	2 µl
Bacterial liquid	2 µl
ddH <sub>2</sub> O	to 20 µl

### Reaction Program


Temperature	Time	Cycles
95°C	3 min	
95°C	15 sec	35
55°C	15 sec	
72°C	15 sec/kb	
72°C	5 min	

**b. Enzyme Digestion Analysis:** According to the experimental design, select the appropriate restriction endonuclease for identification.

**c. Identification of Plasmid Size:** Pick single-cell clones and extract plasmids. Perform electrophoresis to identify the size of plasmids.

**d. Sequencing Analysis:** Directly pick single-cell clones for sequencing identification. Select M13 Forward Primer, M13 Reverse Primer or self-designed primers for sequencing.

## Appendix



M13 Reverse primer site: bases 60 - 76

M13 Forward primer site: bases 209 - 225

Ampicillin resistance ORF (C) : bases 226 - 1,086

pUC origin: bases 1,210 - 1,798

(C) : complementary strand

For complete sequence information of pCE3 Blunt Vector, please visit [www.vazyme.com](http://www.vazyme.com)



**Vazyme**

Tel: +86 25-83772625

Email: [info.biotech@vazyme.com](mailto:info.biotech@vazyme.com)

Web: [www.vazyme.com](http://www.vazyme.com)

Loc: Red Maple Hi-tech Industry Park, Nanjing, PRC