

VAHTS DNA Adapters Set 8 for MGI (Plate)

NMB108

Version 24.1



Product Description

VAHTS DNA Adapters Set 8 for MGI (Plate) is a specialized kit for library preparation based on the MGI sequencers, which is applicable to prepare multi-sample DNA libraries for MGI sequencers. The kit includes 96 different single-index adapters. All the reagents provided in the kit have undergone rigorous quality control and functional testing to ensure the optimal stability and repeatability of library preparation.

Components

Components	NMB108-01 (96 rxns)
DNA Adapter	5 µl each

▲ The concentration of DNA Adapter is 10 µM.

NMB108 Components

	1	2	3	4	5	6	7	8	9	10	11	12
A	01	41	57	65	73	81	89	97	121	25	33	49
B	02	42	58	66	74	82	90	98	122	26	34	50
C	03	43	59	67	75	83	91	99	123	117	35	51
D	04	44	60	68	76	84	92	100	124	28	36	52
E	13	45	61	69	77	85	93	101	125	29	37	53
F	14	46	62	70	78	86	94	102	126	30	38	116
G	15	47	63	71	79	87	95	103	127	114	39	55
H	16	48	64	72	80	88	96	104	128	32	115	56

▲ NMB108-01 layout and combination instructions (It is recommended to use the adapters in the same color box).

Storage

Store at -30 ~ -15°C and ship at ≤0°C.

Applications

It is a specialized kit for VAHTS Universal DNA Library Prep Kit for MGI (Vazyme #NDM607), VAHTS Universal Pro DNA Library Prep Kit for MGI (Vazyme #NDM608), VAHTS Universal Plus DNA Library Prep Kit for MGI (Vazyme #NDM617), VAHTS Universal Plus DNA Library Prep Kit for MGI V2 (Vazyme #NDM627).

Notes

1. Do not pre-mix Adapter, Ligation Buffer and Ligase, as it may cause adapter self-ligation.
2. Do not store this product in an environment higher than room temperature (15 ~ 25°C), otherwise it may lead to a decrease in connection efficiency.

For Research Use Only. Not for use in diagnostic procedures.